

Pan-Succinyl-K Rabbit Polyclonal Antibody For Wb And Elisa Applications

Item Number: CM0000009

Introduction

High-quality rabbit polyclonal antibody targeting pan-succinyl-lysine, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying lysine succinylation and cellular metabolism.

[Learn More](#)

Parameter	Value
Product Name	Pan-Succinyl-K Rabbit pAb
Species	Human
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence containing Succinyl-K.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag