

Pan-Symmetric Di-Methyl Arginine Motif Rabbit Mab For Wb And Elisa Applications

Item Number: CM0000206

Introduction

Pan-Symmetric Di-Methyl Arginine Motif Rabbit monoclonal antibody for WB and ELISA. Reacts with human, mouse, rat, and predicted wide range. Suitable for methylation research.

[Learn More](#)

Parameter	Value
Product Name	Pan-Symmetric Di-Methyl Arginine Motif Rabbit mAb
Species	Human
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence containing Symmetric Di-Methyl R)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag