

Pan-Symmetric Di-Methyl Arginine Rabbit Polyclonal Antibody For Western Blot (Wb) And Elisa

Item Number: CM0000218

Introduction

High-purity pan-Symmetric Di-Methyl Arginine Rabbit Polyclonal Antibody (pAb) for Western Blot (WB) and ELISA, produced against a symmetric dimethylarginine peptide. Offers broad predicted cross-reactivity across human, mouse, and rat for methylation studies.

[Learn More](#)

Parameter	Value
Product Name	pan-Symmetric Di-Methyl Arginine Rabbit pAb
Gene ID (Human)	pan Rme2s
Immunogen	A synthetic peptide corresponding to a sequence containing methylated R.
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag