

Pan-Ubiquitin Remnant Motif (K-E-Gg) Rabbit Polyclonal Antibody For Wb, Elisa

Item Number: CM0000039

Introduction

Anti-pan-Ubiquitin Remnant Motif (K-ε-GG) Rabbit Polyclonal Antibody validated for Western blot and ELISA. Detects ubiquitinated motif across human, mouse, and rat samples. Ideal for ubiquitin research and assay development.

[Learn More](#)

Parameter	Value
Product Name	pan-Ubiquitin Remnant Motif (K-ε-GG) Rabbit pAb
Immunogen	Synthetic peptide containing ubiquitinated Motif (K-ε-GG)
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag