

Par4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Iz0

Item Number: CM0005789



Introduction

High-specificity PAR4 rabbit monoclonal antibody for WB, IHC-P, ELISA. Reacts with human and mouse PAR4, a pro-apoptotic protein for cancer and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	PAR4 Rabbit mAb
Remarks/Alias	PAR4; Par-4
Species	Human
Gene ID (Human)	5074
Gene ID	5074
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 241-340 of human PAR4 (Q96Iz0).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96Iz0
Protein Weight	37kDa
Shipping	Ice bag