

Parg Rabbit Polyclonal Antibody For Wb, Elisa - Q86W56

Item Number: CM0014333



Introduction

PARG Rabbit Polyclonal Antibody validated for Western blotting and ELISA. Recognizes human poly(ADP-ribose) glycohydrolase, a key enzyme in DNA damage repair and ADP-ribosylation metabolism. Predicted molecular weight 111 kDa.

[Learn More](#)

Parameter	Value
Product Name	PARG Rabbit pAb
Remarks/Alias	PARG99; PARG
Species	Human
Gene ID (Human)	8505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 777-976 of human PARG (NP_003622.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86W56
Protein Weight	111kDa
Shipping	Ice bag