

Parp10 Rabbit Polyclonal Antibody For Wb, Elisa - Q53GI7

Item Number: CM0016292



Introduction

Rabbit polyclonal antibody targeting human PARP10 (ARTD10), validated for WB and ELISA. Suitable for studies on DNA damage repair, ADP-ribosylation, and GSK3B signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	PARP10 Rabbit pAb
Remarks/Alias	ARTD10; PARP10
Species	Human
GeneID (human)	84875
GeneID	84875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PARP10 (NP_116178.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53GL7
Protein Weight	110kDa
Shipping	Ice bag