


Parp3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6F1

Item Number: CM0030560



Introduction

Rabbit polyclonal antibody targeting PARP3 (ARTD3) for Western blot, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for studying DNA repair and mono-ADP-ribosylation. Weight 60 kDa.

[Learn More](#)

Parameter	Value
Product Name	PARP3 Rabbit pAb
Remarks/Alias	IRT1; ARTD3; ADPRT3; ADPRTL2; ADPRTL3; PADPRT-3; PARP3
Species	Human
GeneID (Human)	10039
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-533 of human PARP3 (NP_005476.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6F1
Protein Weight	60kDa
Shipping	Ice bag