

Pax3 Rabbit Monoclonal Antibody For Wb, Elisa - P23760

Item Number: CM0007575



Introduction

Rabbit monoclonal anti-PAX3 antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat PAX3, a key transcription factor in neural crest development and myogenesis.

[Learn More](#)

Parameter	Specification
Product Name	PAX3 Rabbit mAb
Remarks/Alias	WS1; WS3; CDHS; HUP2; PAX-3; PAX3
Species	Human
GeneID (Human)	5077
GeneID	5077
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human PAX3 (NP_852122.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23760
Protein Weight	53kDa
Shipping	Ice bag