

Pax4 Rabbit Polyclonal Antibody For Wb, Elisa - O43316

Item Number: CM0026928

Introduction

PAX4 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA. Reacts with human and mouse PAX4, a transcription factor critical for pancreatic beta cell development. Ideal for diabetes and differentiation research.

[Learn More](#)

Parameter	Value
Product Name	PAX4 Rabbit pAb
Remarks/Alias	KPD; MODY9; PAX4
Species	Human
GeneID (Human)	5078
GeneID	5078
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human PAX4 (NP_006184.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43316
Protein Weight	38kDa
Shipping	Ice bag