

Pbxip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Aq6

Item Number: CM0024970



Introduction

Rabbit polyclonal antibody against human PBXIP1 (HPIP), validated for Western blot, IHC-P, and ELISA; detects human, mouse, rat targets; 81 kDa protein; ideal for transcription regulation and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PBXIP1 Rabbit pAb
Remarks/Alias	HPIP; PBXIP1
Species	Human
Gene ID	57326 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human PBXIP1 (NP_065385.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AQ6
Protein Weight	81kDa
Shipping	Ice bag