

Pcgf3 Rabbit Polyclonal Antibody For Wb And Elisa - Q3Knv8

Item Number: CM0025988

Introduction

Rabbit polyclonal antibody against human PCGF3 (RNF3A), a polycomb group protein involved in gene repression and histone modification. Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	PCGF3 Rabbit pAb
Remarks/Alias	RNF3; DONG1; RNF3A; PCGF3
Species	Human
GeneID (Human)	10336
GeneID	10336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human PCGF3 (NP_006306.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q3KNV8
Protein Weight	28kDa
Shipping	Ice bag