

Pcnp Rabbit Pab For Wb, Elisa - Q8Ww12

Item Number: CM0013258

Introduction

Polyclonal antibody targeting human PCNP, a nuclear protein implicated in cell cycle regulation. Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	PCNP Rabbit pAb
Remarks/Alias	PCNP
Species	Human
GeneID (Human)	57092
GeneID	57092
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human PCNP (NP_065090.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WW12
Protein Weight	19kDa
Shipping	Ice bag