

## Pdap1 Rabbit Pab For Wb, If/Icc, Elisa - Q13442

Item Number: CM0028099

### Introduction

Rabbit polyclonal antibody against human PDAP1 (PAP1/HASPP28), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying PDGF signaling and cell growth regulation.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PDAP1 Rabbit pAb                                                                                                         |
| Remarks/Alias    | PAP; PAP1; HASPP28; PDAP1                                                                                                |
| Species          | Human                                                                                                                    |
| GeneID           | 11333                                                                                                                    |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PDAP1 (NP_055706.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human,Mouse,Rat                                                                                                          |
| SWISS            | Q13442                                                                                                                   |
| Protein Weight   | 21kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |