

Pdc Rabbit Pab - P20941

Item Number: CM0017184

Introduction

Phosducin (PDC) rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Recognizes human, mouse, and rat PDC. Recombinant immunogen corresponding to amino acids 1-246 of human PDC. Ideal for visual transduction research.

[Learn More](#)

Parameter	Value
Product Name	PDC Rabbit pAb
Remarks/Alias	PHD; MEKA; PhLP; PhLOP; PDC
Species	Human
GeneID (Human)	5132
GeneID	5132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-246 of human PDC (NP_002588.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P20941
Protein Weight	28kDa
Shipping	Ice bag