

Pdcd2L Rabbit Pab For Wb, Ihc-P, Elisa - Q9Brp1

Item Number: CM0029823

Introduction

PDCD2L Rabbit Polyclonal Antibody for detecting human, mouse, rat PDCD2L. Suitable for Western blot, IHC-P, and ELISA. Ideal for studying ribosome biogenesis and u55 assembly chaperone function.

[Learn More](#)

Parameter	Value
Product Name	PDCD2L Rabbit pAb
Remarks/Alias	PDCD2L
Species	Human
Gene ID	84306
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-358 of human PDCD2L (NP_115722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BRP1
Protein Weight	39kDa
Shipping	Ice bag