

Pde11A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hcr9

Item Number: CM0023679

Introduction

PDE11A Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Shows cross-reactivity with mouse and rat. Suitable for studying cAMP/cGMP hydrolysis and signal transduction in human samples.

[Learn More](#)

Parameter	Value
Product Name	PDE11A Rabbit pAb
Remarks/Alias	PPNAD2; PDE11A
Species	Human
Gene ID (Human)	50940
Gene ID	50940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 712-933 of human PDE11A (NP_058649.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9HCR9
Protein Weight	105kDa
Shipping	Ice bag