


Pde6G Rabbit Pab For Wb, If/Icc, Elisa - P18545

Item Number: CM0019535



Introduction

CamelBio supplies PDE6G Rabbit pAb, a polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat PDE6G, the inhibitory subunit of rod cGMP phosphodiesterase involved in phototransduction.

[Learn More](#)

Parameter	Value
Product Name	PDE6G Rabbit pAb
Remarks/Alias	PDEG; RP57; PDE6G
Species	Human
GeneID (Human)	5148
GeneID	5148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human PDE6G (NP_002593.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18545
Protein Weight	10 kDa
Shipping	Ice bag