

Pdha2 Rabbit Polyclonal Antibody (Wb, If/Icc, Elisa) - P29803

Item Number: CM0027814

Introduction

Rabbit polyclonal antibody targeting human PDHA2 (testis-specific pyruvate dehydrogenase alpha subunit, 43kDa, SWISS P29803). Validated for WB, IF/ICC, ELISA with reactivity in human, mouse, rat. Ideal for studying mitochondrial metabolism and spermatogenesis.

[Learn More](#)

Parameter	Value
Product Name	PDHA2 Rabbit pAb
Remarks/Alias	PDHAL; SPGF70; PDHA2
Species	Human
Gene ID (Human)	5161
Gene ID	5161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 119-388 of human PDHA2 (NP_005381.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29803
Protein Weight	43kDa
Shipping	Ice bag