

Pdhx Rabbit Pab For Wb, Ip, Elisa - 000330

Item Number: CM0027208



Introduction

PDHX Rabbit pAb is a polyclonal antibody validated for WB, IP, and ELISA. It detects human, mouse, and rat PDHX, a mitochondrial protein crucial for pyruvate dehydrogenase complex assembly and energy metabolism.

[Learn More](#)

Parameter	Specification
Product Name	PDHX Rabbit pAb
Remarks/Alias	E3BP; OPDX; PDX1; proX; DLDBP; PDHXd; PDHX
Species	Human
GeneID (Human)	8050
GeneID	8050
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00330
Protein Weight	54kDa
Shipping	Ice bag