

Pdk4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16654

Item Number: CM0028757



Introduction

Anti-PDK4 rabbit polyclonal antibody validated for Western blotting, IHC-P, and ELISA. Detects human, mouse, and rat PDK4, a mitochondrial kinase regulating glucose and fatty acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	PDK4 Rabbit pAb
Remarks/Alias	PDK4
Species	Human
GeneID (Human)	5166
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 312-411 of human PDK4 (NP_002603.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16654
Protein Weight	46kDa
Shipping	Ice bag