

Pdk4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16654

Item Number: CM0016813



Introduction

High-quality PDK4 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA applications. Reacts with human, mouse, rat PDK4 (46kDa). Ideal for studying glucose and fatty acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	PDK4 Rabbit pAb
Remarks/Alias	PDK4
Species	Human
GeneID (human)	5166
GeneID	5166
Immunogen	Synthetic peptide corresponding to amino acids 102-411 of human PDK4 (NP_002603.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16654
Protein Weight	46kDa
Shipping	Ice bag