

Pdxk Rabbit Mab For Wb, If/Icc, Elisa - O00764

Item Number: CM0011704



Introduction

PDXK rabbit monoclonal antibody validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for studying vitamin B6 metabolism, PLP-dependent enzymes, and Charcot-Marie-Tooth disease type 6C.

[Learn More](#)

Parameter	Value
Product Name	PDXK Rabbit mAb
Remarks/Alias	PKH, PNK, HMSN6C, PRED79, C21orf97, HEL-S-1a, C21orf124
Species	Human
GeneID (human)	8566
GeneID	8566
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 213-312 of human PDXK (NP_003672.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O00764
Protein Weight	27kDa/32kDa/35kDa
Shipping	Ice bag