

Pe Anti-Cd94 Monoclonal Antibody For Flow Cytometry - O54707

Item Number: CM0009437

Introduction

PE-conjugated rabbit monoclonal antibody against mouse CD94 (Klrp1), validated for flow cytometry. Targets NK cell inhibitory receptor for immune checkpoint and immuno-oncology research. Cross-reacts with mouse. 21kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD94 mAb
Remarks/Alias	CD94
Species	Mouse
GeneID (Human)	3824
GeneID	16643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-179 of mouse CD94 (NP_034784.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O54707
Protein Weight	21kDa
Shipping	Ice bag