

Pe Anti-Human Cd15/Ssea-1 Monoclonal Antibody For Flow Cytometry - P22083

Item Number: CM0004915

Introduction

This PE-conjugated mouse monoclonal antibody specifically targets human CD15/SSEA-1 (FUT4), a glycosyltransferase that generates Lewis X antigen. Ideal for flow cytometry (FC) detection. Swiss accession: P22083.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD15/SSEA-1 mAb
Remarks/Alias	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
Species	Human
GeneID (human)	2526
GeneID	2526
Immunogen	Synthetic peptide
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22083
Protein Weight	59kDa
Shipping	Ice bag