

Pe Anti-Human Cd171/L1Cam Monoclonal Antibody For Flow Cytometry - P32004

Item Number: CM0011359



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD171/L1CAM, validated for flow cytometry. Supports research on neural cell adhesion, brain development, and L1CAM-associated neurological disorders.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD171/L1CAM mAb
Remarks/Alias	S10; HSAS; HYCX; MASA; MIC5; SPG1; CAML1; CD171; HSAS1; N-CAML1; NCAM-L1; N-CAM-L1
Species	Human
GeneID (Human)	3897
GeneID	3897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-1120 of human CD171/L1CAM (NP_000416.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32004
Protein Weight	140kDa/139kDa
Shipping	Ice bag