


Pe Anti-Human Cd209/Dc-Sign Monoclonal Antibody For Flow Cytometry - Q9Nnx6

Item Number: CM0009717



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD209/DC-SIGN, validated for flow cytometry. Suitable for studying dendritic cell function, pathogen recognition, and immune response in human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD209/DC-SIGN mAb
Remarks/Alias	CDSIGN; CLEC4L; DC-SIGN; DC-SIGN1; hDC-SIGN
Species	Human
Gene ID (human)	30835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 59-404 of human CD209/DC-SIGN (NP_066978.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NNX6
Protein Weight	4kDa/18kDa/30kDa/33kDa/35kDa/37kDa/41kDa/43kDa/45kDa
Shipping	Ice bag