


Pe Anti-Human Cd22 Monoclonal Antibody For Flow Cytometry - P20273

Item Number: CM0004173



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD22 (Siglec-2), validated for flow cytometry. Ideal for B-cell immunophenotyping, lymphoma research, and B-cell receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD22 mAb
Remarks/Alias	SIGLEC2; SIGLEC-2
Species	Human
GeneID (Human)	933
GeneID	933
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20273
Protein Weight	95kDa
Shipping	Ice bag