

Pe Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry - P26842

Item Number: CM0006181



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD27, a costimulatory immune-checkpoint receptor on T-cells, B-cells, and NK-cells. Validated for flow cytometry to study CD27-mediated T-cell activation and memory formation.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
Gene ID (Human)	939
Gene ID	939
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag