

Pe Anti-Human Cr2/Cd21 Monoclonal Antibody For Fc - P20023

Item Number: CM0004746



Introduction

PE-conjugated rabbit monoclonal antibody against human CR2/CD21, validated for flow cytometry. Targets complement receptor type 2 (C3d/EBV receptor), a B-cell coreceptor and viral entry mediator. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CR2/CD21 mAb
Remarks/Alias	CR; C3DR; CD21; CVID7; SLEB9
Species	Human
Gene ID (Human)	1380
Gene ID	1380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-971 of human CD21 (NP_001868.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20023
Protein Weight	110kDa/113kDa/118kDa/119kDa
Shipping	Ice bag