

Pe Anti-Met/C-Met Monoclonal Antibody For Flow Cytometry - P08581

Item Number: CM0010444



Introduction

This PE-conjugated rabbit monoclonal antibody targets human Met/c-Met, validated for intracellular flow cytometry (FC). Recombinant immunogen (aa 25-932) ensures specificity, supporting studies in oncology, development, and receptor signaling.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human Met/c-Met mAb
Remarks/Alias	DA11, HGFR, AUTS9, RCCP2, c-Met, DFNB97
Species	Human
Gene ID (Human)	4233
Gene ID	4233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-932 of human Met/c-Met (NP_000236.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P08581
Protein Weight	156kDa
Shipping	Ice bag