

Pe Anti-Mouse Cd11C Mab - Q9Qxh4

Item Number: CM0005831

Introduction

PE-conjugated Armenian hamster monoclonal antibody targeting mouse CD11c (Integrin alpha-X). Designed for flow cytometry applications to study inflammatory responses, monocyte adhesion, and erythrophagocytosis in mouse models.

[Learn More](#)

Parameter	Value
Product Name	PE anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag