

Pe Anti-Mouse Cd69 Monoclonal Antibody For Flow Cytometry - P37217

Item Number: CM0008528

Introduction

Rabbit monoclonal antibody conjugated to PE, specific for mouse CD69 (SWISS P37217). Validated for flow cytometry. Ideal for studying T-cell activation and immune cell homeostasis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
Gene ID (Human)	969
Gene ID	12515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-199 of mouse CD69 (NP_001028294.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag