


Pe Anti-Mouse Mhc Class Ii (I-A/I-E) Monoclonal Antibody For Flow Cytometry - P18468 / P06342 / Q3U060 / P14435

Item Number: CM0009858



Introduction

PE-conjugated rat monoclonal antibody targeting mouse MHC Class II I-A/I-E molecules. Validated for flow cytometry, suitable for immunology and immuno-oncology research. Specificity for P18468, P06342, Q3U060, P14435.

[Learn More](#)

Parameter	Value
Product Name	PE Rat anti-Mouse MHC Class II (I-A/I-E) mAb
Remarks/Alias	Eb; Ia4; H2Eb; Ia-4; H-2Eb/IAb; Ia2; Ia-2; Abeta; H-2Ab; H2-Ab; Rmcs1; I-Abeta; Ia5; Ia-5; H-2Eb2; Iebeta2; A130038H09Rik; Ia1; H2Aa; Ia-1; H-2Aa; Aalpha; IAalpha; I-Aalpha
Species	Mouse
GeneID	14969/14961/381091/14960
Immunogen	Recombinant protein of mouse MHC Class II (I-A/I-E)
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18468/P06342/Q3U060/P14435
Protein Weight	30kDa/32kDa/26kDa
Shipping	Ice bag