

Pe Anti-Mouse Tcr γ/δ Monoclonal Antibody For Flow Cytometry - Tcrdelta

Item Number: CM0007817

Introduction

PE-conjugated Armenian hamster monoclonal antibody against mouse TCR γ/δ . Validated for flow cytometry, suitable for studying $\gamma\delta$ T cell populations in mouse immune system research.

[Learn More](#)

Parameter	Specification
Product Name	PE anti-Mouse TCR γ/δ mAb
Remarks / Alias	Tcrdelta; TCRGV1S1; TCRGV2S1; TCRGV3S1; TCRGV5S3
Species	Mouse
Gene ID	110066/110067
Immunogen	Recombinant protein of mouse TCR γ/δ
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
Shipping	Ice bag