

Pe Anti-Mouse/Human Cd55 Monoclonal Antibody For Flow Cytometry - Q61475

Item Number: CM0009296

Introduction

PE-conjugated rabbit monoclonal antibody against mouse/human CD55 (DAF) validated for flow cytometry. Targets the complement decay-accelerating factor, a GPI-anchored protein that inhibits C3/C5 convertases. Ideal for immuno-oncology and complement research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse/Human CD55 mAb
Remarks/Alias	Daf; Daf1; Daf-GPI; GPI-DAF
Species	Mouse
GeneID (Human)	1604
GeneID	13136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-361 of mouse CD55/DAF (NP_034146.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	Q61475
Protein Weight	42KDa
Shipping	Ice bag