


Pe Anti-Vegfr-3/Flt-4 Rabbit Monoclonal Antibody For Flow Cytometry - P35916

Item Number: CM0011762



Introduction

PE-conjugated rabbit monoclonal antibody targeting human VEGFR-3/FLT-4. Recombinant immunogen (aa 25-776). Validated for flow cytometry detection in human samples. Ideal for lymphangiogenesis and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human VEGFR-3/FLT-4 mAb
Remarks/Alias	PCL, CHTD7, FLT-4, FLT41, LMPH1A, LMPHM1, VEGFR3, VEGFR-3
Species	Human
GeneID (Human)	2324
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-776 of human VEGFR3/FLT4 (NP_002011.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P35916
Protein Weight	153 kDa
Shipping	Ice bag