

Pe-Cf594 Rabbit Anti-Human Bcl6 Mab - P41182

Item Number: CM0009765



Introduction

Conjugated monoclonal antibody targeting human BCL6 for intracellular flow cytometry. Recognizes BCL6, a transcriptional repressor essential for germinal center reactions, with applications in immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	PE-CF594 Rabbit anti-Human BCL6 mAb
Remarks/Alias	BCL5; LAZ3; BCL6A; ZNF51; ZBTB27
Species	Human
GeneID (Human)	604
GeneID	604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-350 of human BCL6 (NP_001124317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P41182
Protein Weight	79kDa
Shipping	Ice bag