

Pe-Cf594 Rabbit Anti-Human/Monkey Cd69 Mab For Flow Cytometry - Q07108

Item Number: CM0008663



Introduction

PE-CF594-conjugated rabbit monoclonal antibody targeting human and monkey CD69, validated for flow cytometry. Detects early activation antigen on T-cells, B-cells, NK cells; 23kDa; cross-reacts with cynomolgus.

[Learn More](#)

Parameter	Value
Product Name	PE-CF594 Rabbit anti-Human/Monkey CD69 mAb
Remarks/Alias	AIM; EA1; MLR-3; CLEC2C; GP32/28; BL-AC/P26
Species	Human
GeneID (human)	969
GeneID	969
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 78-199 of human CD69 (NP_001772.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	Q07108
Protein Weight	23kDa
Shipping	Ice bag