

Pe-Cf594 Rabbit Anti-Human/Monkey Hla-Dr Monoclonal Antibody For Flow Cytometry - P01903

Item Number: CM0009779



Introduction

PE-CF594-conjugated rabbit monoclonal antibody against human and cynomolgus monkey HLA-DR alpha chain (P01903), optimized for flow cytometry (FC). Enables detection of MHC class II-expressing antigen-presenting cells in immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	PE-CF594 Rabbit anti-Human/Monkey HLA-DR mAb
Remarks/Alias	HLA-DRA1
Species	Human
GeneID (Human)	3122
GeneID	3122
Immunogen	Recombinant protein of human HLA-DR.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Cynomolgus
SWISS	P01903
Protein Weight	29kDa
Shipping	Ice bag