

Pe-Cf594 Rabbit Anti-Human/Mouse Pdgfra/Cd140A Monoclonal Antibody For Fc - P16234

Item Number: CM0005042



Introduction

PE-CF594-conjugated rabbit monoclonal antibody against human and mouse PDGFR α /CD140a, validated for flow cytometry. Suitable for studying PDGF receptor signaling, cell proliferation, embryonic development, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PE-CF594 Rabbit anti-Human/Mouse PDGFR α /CD140a mAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2
Species	Human/Mouse
Gene ID (Human)	5156
Gene ID	5156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-524 of human/mouse PDGFR α (NP_006197.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human, Mouse
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag