

Pe-Conjugated Anti-Mouse Cd22 Monoclonal Antibody For Flow Cytometry - P35329

Item Number: CM0009149

Introduction

This PE-conjugated rabbit anti-mouse CD22 monoclonal antibody is designed for flow cytometry. It binds the B-cell receptor CD22 (Siglec-2) and facilitates B-cell identification and functional analysis in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
GenelD (human)	933
GenelD	12483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa/93kDa/97kDa
Shipping	Ice bag