

Pe Mouse Anti-Human CD11B Monoclonal Antibody For Fc - P11215

Item Number: CM0004505



Introduction

PE-conjugated mouse monoclonal antibody targeting human CD11b (Integrin alpha-M, MAC-1) for flow cytometry. Specific for human CD11b, a 127 kDa leukocyte adhesion receptor involved in innate immunity and inflammation.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD11b mAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
GeneID (Human)	3684
GeneID	3684
Immunogen	A synthesized peptide derived from human CD11b
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag