

Pe Mouse Anti-Human Cd11C Monoclonal Antibody For Flow Cytometry - P20702

Item Number: CM0003993



Introduction

PE-conjugated mouse monoclonal antibody specifically targeting human CD11c (Integrin alpha-X, ITGAX). Validated for flow cytometry (FC) and cross-reactive with human samples. Useful for studying monocyte adhesion, chemotaxis, and erythrophagocytosis.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD11c mAb
Remarks/Alias	CD11C; SLEB6
Species	Human
Gene ID (Human)	3687
Immunogen	Recombinant protein of human CD11c.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20702
Protein Weight	128kDa
Shipping	Ice bag