

## Pe Mouse Anti-Human Cd15/Ssea-1 Monoclonal Antibody For Fc - P22083

Item Number: CM0012294

### Introduction

PE-conjugated mouse monoclonal antibody targeting human CD15/SSEA-1 (FUT4) epitope, validated for flow cytometry. Suitable for research in leukocyte biology, stem cell characterization, and selectin-mediated adhesion studies.

[Learn More](#)

| Parameter        | Value                                            |
|------------------|--------------------------------------------------|
| Product Name     | PE Mouse anti-Human CD15/SSEA-1 mAb              |
| Remarks/Alias    | LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV   |
| Species          | Human                                            |
| GeneID (Human)   | 2526                                             |
| GeneID           | 2526                                             |
| Immunogen        | Synthetic peptide derived from human CD15/SSEA-1 |
| Source           | Mouse                                            |
| Category         | Monoclonal Antibodies                            |
| Application      | FC                                               |
| Cross Reactivity | Human                                            |
| SWISS            | P22083                                           |
| Protein Weight   | 59kDa                                            |
| Shipping         | Ice bag                                          |