

Pe Mouse Anti-Human Cd1A Monoclonal Antibody For Fc - P06126

Item Number: CM0005154



Introduction

PE-conjugated mouse monoclonal antibody targeting human CD1a, a lipid antigen-presenting molecule, validated for flow cytometry. Ideal for studying NKT cell activation and dendritic cell function in cancer and infectious diseases.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD1a mAb
Remarks/Alias	R4; T6; CD1; FCB6; HTA1
Species	Human
GeneID (Human)	909
GeneID	909
Immunogen	Synthetic peptide A synthesized peptide derived from human CD1a.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06126
Protein Weight	37kDa
Shipping	Ice bag