

Pe Mouse Anti-Human Cd3 Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0011580



Introduction

PE-conjugated mouse monoclonal antibody targeting human CD3 epsilon chain (UniProt: P07766). Suitable for flow cytometry-based T-cell detection and immunophenotyping. This reagent enables identification of T lymphocytes in various sample types. Species reactivity: Human.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag