

Pe Mouse Anti-Human Cd42A Monoclonal Antibody For Fc - P14770

Item Number: CM0010870



Introduction

PE-conjugated mouse monoclonal antibody targeting human CD42a (GPIX). For flow cytometry (FC) analysis of platelet glycoprotein IX, a subunit of the GPIb-V-IX von Willebrand factor receptor. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD42a mAb
Remarks/Alias	GPIX; CD42a
Species	Human
Gene ID	2815
Immunogen	Recombinant protein of human CD42a.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS Accession	P14770
Protein Weight	19kDa
Shipping	Ice bag