

Pe Mouse Anti-Human Cd45Ro Monoclonal Antibody For Flow Cytometry - P08575

Item Number: CM0012215



Introduction

PE-conjugated mouse monoclonal antibody targeting human CD45RO, a receptor-type tyrosine-protein phosphatase C involved in T-cell activation. Validated for flow cytometry (FC) applications, recognizing the human CD45RO isoform.

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD45RO mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
Immunogen	Synthetic peptide
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575
Protein Weight	136kDa
Shipping	Ice bag