

Pe Mouse Anti-Human Cd56 Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0011175



Introduction

PE-conjugated mouse monoclonal antibody against human CD56 (NCAM1), validated for flow cytometry (FC). A reliable tool for studying natural killer cells, neural adhesion, and viral entry mechanisms. - P13591

[Learn More](#)

Parameter	Value
Product Name	PE Mouse anti-Human CD56 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
Gene ID	4684
Immunogen	Synthetic peptide A synthesized peptide derived from human CD56
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	40kDa/73kDa/80kDa/83kDa/93kDa/95kDa
Shipping	Ice bag